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
The patient is markedly hypoalbuminemic. Severe hypoalbuminemia may be associated with decreased production (i.e. compromised hepatic function) and/or increased loss (e.g., PLE, PLN). PLN or glomerulopathy is the most commonly recognized cause of marked hypoalbuminemia (nephrotic syndrome). Albumin concentrations of 1.5 g/dl or less are generally associated with dependent edema and/or ascites. The alpha-2 globulin fraction is elevated. The increase is non-specific, as this band contains a majority of the protease inhibitors associated with

acute inflammation. An elevation will generally accompany any inflammatory process.

The gamma globulin fraction is elevated, characterized by a **broad-based polyclonal band, resulting from a mixture of immunoglobulins associated with an immune response**. Potential causes include suppurative disease, chronic infectious disease (bacterial, viral, protozoal, rickettsial, fungal etc), immune mediated disease, chronic granulomatous disease, or inflammation secondary to neoplasia, but not secretory round cell neoplasia (e.g. multiple myeloma or B cell lymphoproliferative disease), as there is no evidence of a monoclonal gammopathy. Correlate with clinical findings, and other diagnostic tests. **FIP is a common cause of this pattern in young cats; however, consider also other chronic infectious agents, septic processes or inflammation involving viscera.**

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